

PRODUCT INFORMATION

Model: RLDF0205 Bench top/Wall mounted Refrigerator

Type: Labcold Pharmacy/Vaccine Refrigerator

Designed for bench top storage, this professional pharmacy fridge is ideal for the surgery, pharmacy or ward. The integral lock ensures that contents are kept secure and the audio alarm makes sure that they are kept safely. Completely complies with all current guidance regarding the safe storage of vaccines. Can be wall mounted on integral brackets if required. This compact fridge auto defrosts for convenience and like the rest of the Labcold range is CFC, HCFC and ammonia free.



Capacity: 66 litres

Construction: White enamel exterior, metal interior

Shelves Supplied: 3 wire mesh

Exterior Dimensions (mm): 720H x 445W x 475D

Features:

Integrated digital display of temperature with max/min recording

Audible alarm with mute facility

Fan assisted air circulation

Integral wall mounting brackets

Door lock

CFC and HCFC free and Ammonia free

Accurate temperature control: +2°C to +8°C

Automatic defrost and removal of condensate

Order Reference: RLDF0205