

PRODUCT INFORMATION

Model: RLDF0505

Type: Labcold Pharmacy Refrigerator

A purpose-built refrigerator for laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with NEW UK Department of Health "Green Book" requirements, plus recommendations in Royal Pharmaceutical Society of Great Britain's "Code of Ethics and Standards".



Capacity: 150 litres

Construction: White enamel exterior, white painted steel interior.

Shelves supplied: 3, plastic coated wire mesh.

Exterior dimensions (mm): 820H x 600W x 600D

Features:

Integral digital temperature display

Max/min temperature record with reset

Audible alarm with mute facility

Fan assisted air circulation

Internal light

Door lock, with two keys

Integral door handle

CFC, HCFC and Ammonia free

Accurate temperature control: +2°C to +8°C

Automatic defrost, with condensate removal

Order Reference: RLDF 0505

Optional Accessories: Pharmacy Drawer System,

Order Reference: PHARMDRAWER01