

## PRODUCT INFORMATION

Model: RLDF1005, RLDG1005 (Glass Door)

Type: Labcold Pharmacy/Vaccine Refrigerator

**A purpose-built refrigerator for laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with NEW UK Department of Health “Green Book” requirements, plus recommendations in Royal Pharmaceutical Society of Great Britain’s “Code of Ethics and Standards”.**



**Capacity: 300 litres**

**Construction: White enamel exterior, white painted steel interior**

**Shelves Supplied: 6 wire mesh**

**Exterior Dimensions (mm): 1487H x 600W x 650D**

**Features:**

**Digital display of temperature with max/min record**

**Audible alarm with mute facility**

**Max/min temperature memory and display**

**Fan assisted air circulation**

**Internal light**

**Door lock**

**CFC and HCFC free and Ammonia free**

**Accurate temperature control: +2°C to +8°C**

**Integrated digital temperature display and controller**

**Automatic Defrost**

**Order Reference: RLDF 1005**

**RLDG 1005**

**Optional Accessories: Pharmacy Drawer System,**

**Order Reference: PHARMDRAWER02**

**Labcold Limited, The Loddon Centre, Wade Road, Basingstoke, RG24 8FL, UK**