

## PRODUCT INFORMATION

**Model: RLDG0505**

**Type: Labcold Pharmacy Refrigerator (Glass Door)**

A purpose-built refrigerator for laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with NEW UK Department of Health "Green Book" requirements, plus recommendations in Royal Pharmaceutical Society of Great Britain's "Code of Ethics and Standards".



**Capacity: 150 litres**

**Construction: White enamel exterior, white painted steel interior.**

**Shelves supplied: 3, plastic coated wire mesh.**

**Exterior dimensions (mm): 820H x 600W x 600D**

**Features:**

- Integral digital temperature display**
- Max/min temperature record with reset**
- Audible alarm with mute facility**
- Fan assisted air circulation**
- Internal light**
- Door lock, with two keys**
- Integral door handle**
- CFC, HCFC and Ammonia free**
- Accurate temperature control: +2°C to +8°C**
- Automatic defrost, with condensate removal**

**Order Reference: RLDG 0505**

**Optional Accessories: Pharmacy Drawer System,**

**Order Reference: PHARMDRAWER01**