

## **PRODUCT INFORMATION**

**Model: RLFR1804**

**Type: Labcold Fan Circulated Basic Refrigerator**

For laboratory and general use where large capacity combined with value for money is required and there is no risk that volatile, flammable materials or components might be stored.



**Capacity: 505 litres (gross)**

**Construction: White enamelled steel exterior and interior**

**Shelves supplied: 4; Plastic coated wire mesh**

**Exterior dimensions (mm): 1710 high x 765 wide x 725 deep**

**Features:**

**Fan circulated interior**

**Door lock**

**CFC and HCFC free**

**Temperature set-point range: +2 to +5°C**

**Analogue temperature display**

**Automatic defrost**

**Order reference: RLFR 1804**