

PRODUCT INFORMATION

Model: RLHF13043

Type: Labcold Fan Circulated Laboratory Refrigerator

For laboratory and general use where improved temperature uniformity is required and there is no risk that volatile, flammable materials or components might be stored.



Capacity: 390 litres (gross)

Construction: White coated steel, fully moulded plastic interior.

Shelves Supplied: 6; Plastic coated wire mesh

Exterior Dimensions (mm): 1874 high x 600 wide x 615 deep

Features:

Fan circulated interior

Door Lock

Easy-grip door handle

CFC and HCFC free

Temperature set-point range: +2°C to +8°C

Temperature display

Automatic Defrost

Order Reference: RLHF13043