

## **PRODUCT INFORMATION**

**Model: RLPR05042A**

**Type: Labcold Sparkfree Laboratory Refrigerator**

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



**Capacity: 160 Litres (gross)**

**Construction: White enamel exterior, fully moulded plastic interior**

**Shelves Supplied: 3 x toughened glass shelves, 1 x drawer.**

**Exterior Dimensions (mm): 850H x 550W x 620D**

**Features:**

- Sparkfree Interior**
- Removable door pockets**
- Door lock**
- Easy-grip door handle**
- CFC and HCFC free**
- Temperature set-point range: 0 to +10°C**
- Low energy operation**
- Automatic defrost**

**Order Reference: RLPR 05042A**