

PRODUCT INFORMATION

Model: RLPR09043A

Type: Labcold Sparkfree Laboratory Refrigerator

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



Image for illustrative purposes only

Capacity: 285 litres (gross)

Construction: White enamel exterior, fully moulded plastic interior

Shelves supplied: 5, glass

Exterior dimensions (mm): 1435 high x 600 wide x 625 deep

Features:

Sparkfree interior

Door lock

CFC and HCFC free

Temperature set-point range: +2 to +8°C

Automatic defrost

Order reference: RLPR 09043A