

PRODUCT INFORMATION

Model: RLVF02203

Type: Labcold Sparkfree Laboratory Freezer (-20°C)

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



Capacity: 36 Litres (gross)

Construction: White enamel exterior, fully moulded plastic interior.

Shelves supplied: 1, wire mesh

Exterior dimensions (mm): 515 high x 440 wide x 500 deep

Features:

Sparkfree interior

Counter-top use

CFC and HCFC free

Temperature set-point range: -10 to -26°C

Operator controlled defrost

Order reference: RLVF 02203