

PRODUCT INFORMATION

Model: RLVF04202A

Type: Labcold Sparkfree Laboratory Freezer (-20°C)

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



Capacity: 110 Litres (gross)

Construction: White enamel exterior, fully moulded plastic interior.

Internal arrangements: Four plastic drawers separated by fixed refrigerated shelves.

Exterior Dimensions (mm): 850 high x 545 wide x 565 deep

Features:

- Sparkfree Interior**
- Removable freezer drawers**
- CFC and HFC free**
- Temperature set-point range: -18 to -23°C**
- Low energy operation**
- Operator controlled defrost**

This model is an alternate substitute product for model RLVF04202

Order Reference: RLVF 04202A