

## **PRODUCT INFORMATION**

**Model: RLVF14201**

**Type: Labcold Sparkfree Laboratory Freezer (-20°C)**

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



**Capacity: 375 litres (gross)**

**Construction: White enamel exterior, fully moulded plastic interior.**

**Shelves supplied: 6; wire mesh, fixed position**

**Exterior dimensions (mm): 1850 high x 595 wide x 640 deep**

**Features:**

**Sparkfree interior**

**Door lock**

**Easy-grip door handle**

**CFC and HCFC free**

**Temperature set-point range: -18 to -23°C**

**Digital temperature display**

**Operator controlled defrost**

**Order reference: RLVF 14201**

**Labcold Limited, The Loddon Centre, Wade Road, Basingstoke, RG24 8FL, UK**