

PRODUCT INFORMATION

Model: RLVF1825

Type: Labcold Large Capacity Basic Freezer

For laboratory and general use where large capacity combined with value for money is required and there is no risk that volatile, flammable materials or components might be stored. Also suitable for storage of certain temperature sensitive materials such as enzymes.



Capacity: 505 litres (gross)

Construction: White enamelled steel exterior and interior.

Shelves supplied: 6; fixed position, refrigerated

Exterior dimensions (mm): 1710 high x 765 wide x 725 deep

Features:

Door lock

CFC and HCFC free

Temperature set-point range: -18 to -25°C

Analogue temperature display

High temperature warning light

Operator controlled defrost

Order reference: RLVF 1825