

PRODUCT INFORMATION

Model: RLVL03203

Type: Labcold Basic Freezer

For laboratory and general use where reasonable performance and low cost are important.



Capacity: 70 Litres (gross)

Construction: White enamel exterior, fully moulded plastic interior.

Internal arrangement: Three drawers sections separated by fixed refrigerated shelves.

Exterior Dimensions (mm): 845 high x 494 wide x 494 deep

Features:

Two removable freezer drawers

CFC and HCFC free

Temperature set-point range: -15 to -20°C

Operator controlled defrost

Order Reference: RLVL 03203