

PRODUCT INFORMATION

Model: RPFR05043 (Solid Door)

RPFG05043 (Glass Door)

Type: Labcold Advanced Pharmacy Refrigerator

For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.



Capacity: 150 litres (gross)

Construction: AISI 304 grade stainless steel interior and exterior.

Shelves Supplied: 2, Plastic coated wire mesh 530 x 325mm.

Exterior Dimensions (mm): 820 high x 600 wide x 620 deep.

Features:

Microprocessor temperature control with digital display to 1°C

Hi/Lo temperature alarm - audible

Fan circulated interior

Door Lock

CFC and HFC free

Temperature set-point range: +2 to +8°C

Automatic defrost with condensate removal

Maximum ambient temperature: +43°C

Order Reference: RPFR 05043

RPFG 05043